



## 4200 Series Pneumatic Stools

**Tradition and Innovation** – Since 1858, visionary thinking and a fascination with technology have guided us to develop innovative products of outstanding reliability: Anticipating trends to improve the quality of life.

 **RELIANCE  
INTERNATIONAL**

# 4200 Series Pneumatic Stools

## Finished for Durability

- Select from 14 easy-to-maintain vinyl colors. Other vinyl colors available by request.
- All exposed metal, with exception of the aluminum base, is chrome-plated steel.
- Base and legs finished in tough, black epoxy enamel.
- Hooded dual-wheel casters.
- 5 legs for stability.

## Optional Foot Ring

Order a foot ring with your stool, or at a later date.

## Tech Specs

Base diameter: 22.5"



	4240	4246
Seat adjustment	Pneumatic Assist	Pneumatic Assist
Seat diameter	15"	15"
Seat construction	2" polyfoam	2" polyfoam
Seat swivels	360°	360°
Five legs for extra stability	Yes	Yes
Height range	19.50-27"	19.50-27"
Shipping weight	23 lbs	26 lbs