



5300 Series Hydraulic Stools

Tradition and Innovation – Since 1858, visionary thinking and a fascination with technology have guided us to develop innovative products of outstanding reliability: Anticipating trends to improve the quality of life.

 **RELIANCE
INTERNATIONAL**

5300 Series Hydraulic Stools

Exceptional Stability

Finish and Upholstery

- Select from 14 easy-to-maintain vinyl colors. Other vinyl colors available by request.
- All exposed metal, with exception of the aluminum base, is chrome-plated steel.
- Rugged 2½" hard rubber casters.

Easily Adjustable Height Mechanism

Lift the control ring located under the seat. The stool's pneumatic device will help lift the seat as you ease your weight off the stool. Vertical travel is 6.25". Seat height ranges from 18.50" to 24.75". Footrest ring available as an option.

Tech Specs

Base diameter 19¾". Removable footrest: 16" diameter.



	5340	5346	5348	5356
Seat dimensions	15" circular	15" circular	16 x 14" sloped waterfall construction supports thighs	15" circular
Seat construction	2" polyfoam	2" polyfoam	2" polyfoam	2" polyfoam
Swivel	360°	360°	360°	360°
Five legs for extra stability	Yes	Yes	Yes	Yes
Circular control ring for pneumatic lift	Yes	Yes	Yes	Yes
Upholstered 3-way, wrap around arm	No	No	No	Yes, exceptional abdominal, back or arm support
Shipping Weight	26 lbs	32 lbs	32 lbs	47 lbs