



7800 Ophthalmic Instrument Stand

Tradition and Innovation – Since 1858, visionary thinking and a fascination with technology have guided us to develop innovative products of outstanding reliability: Anticipating trends to improve the quality of life.

7800

Ophthalmic Instrument Stand

The Perfect Match for all Reliance Chairs

The 7800 IC (with instrument charging console) and the 7800 NC (no console) are designed for compatibility with any Reliance examination chair. Its sleek, contoured design offers efficient operation and easy access to patients.

Counterbalanced Slit Lamp Arm

- Lock and adjust with a single hand. Arm travels 12" and counterbalances for instruments 12-60 lbs.
- A slit lamp arm-mounted switch raises and lowers the chair.

Upper Refractor Arm

10" range of movement; delivers instruments up to 20 lbs.

Optional Third Instrument Arm

- 20- 50 lb. capacity; 10" of vertical travel.
- Compatible with 5399 Monitor Arm.

Instrument Wells (IC only)

Console is equipped with three rechargeable instrument wells that hold hand-held instruments. It also features raise, lower and auto down switches for the examination chair.

Backlit Membrane Switches (IC only)

Membranes last 10-times longer than mechanical switches and increase safety with lower voltages. They conform the chart projector, overhead lamp and binding post for an indirect ophthalmoscope, and raise and lower the exam chair.

Halogen Lamp

Long-lasting bulb provides bright, direct light for minor procedures. Arm adjusts without tools. Low voltage for safety.

CSA Mark of Quality

Indicates that the 7800 has passed all rigorous CSA,C/US tests, affirming its safety and durability and achieving certificate of compliance.

OPTIONS

Wheelchair Access

Lower slit lamp arm facilitates exams for wheelchair bound patients. Cost eligible for ADA tax credit.

Eclipse Room Light Control

The Eclipse activates all room lighting from the stand, or a small console. Lighting adjusts based on instrument moves or quick manual adjustments.

Tech Specs

Height to console: 36" Height to pole top: 71" Base: 21 x 21" Shipping weight: 292 lbs. Standard electrical: 115 volts, 60 Hz, 3 amperes, optional 230 volts, 50 Hz, 1.5 amperes.

