



GOE Digital PD Meter

Now with LED light-No more changing bulbs!

Precision instrument used to measure the distance between a patients pupils for properly fitting eye glass lenses. The Optician or Optometrist simply aligns the instruments sliders with the reflection of the LED target light on the patients cornea for measurements accurate within .5mm!

Features include:

- Measurement distance range from 35cm to Infinity
- PD Measurement range, Monocular (one eye) 24.5– 38mm, Binocular (Both eye) 49mm– 76mm
- 3 LCD Displays for: Right Monocular, Left Monocular and combined Binocular measurements
- Measurements are accurate to within .25mm per eye
- Automatically powers OFF when not in use
- 10.25" long x 6.5" wide x 2.5" tall
- Instrument weight is only 2 lbs with a 3 lbs. Shipping weight
- LED Bulb: Brighter target light and never needs changing!
- Power requirements: 3 AAA Batteries-Included
- GOE Stock # GOE/414N8L
- CE Approved



Full 1 Year Limited
Warranty

List Price \$425.00	Special Price \$295.00
-------------------------------	----------------------------------



Grimes Optical Equipment Company

Your Best Finishing Equipment Value!

1201 W. 6th Ave.
Amarillo, Texas 79101
U.S.A.

Phone: 800-749-8427
Phone: 806-342-4200
Fax: 806-342-4201
Email: info@grimesoptical.com